Coming Home®

ShowerAid


Skin Barrier/
Adhesive Remover

Skin Barrier/
Adhesive Remover

Figure 1

1.1

1.2

1.3

Figure 2

2.1

2.2

2.3




Figure 3
StabiliBase™, StabiliBase™ OptiDerm™


3.1


3.2


3.3


3.4


3.5

3.6

3.7

3.8




Figure 4

FlexiDerm™, XtraBase®, OptiDerm™


4.1

4.2

4.3

4.4



Figure 5


5.1


5.2


5.3

5.4


5.5


5.6


5.7

5.8

5.9

Disclaimer

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Patents and trademarks

Provox® and Coming Home® are registered trademarks owned by Atos Medical AB, Sweden.

Micron HME™, FlexiDerm™, OptiDerm™, StabiliBase™, XtraHME™, XtraMoist™, XtraFlow™ and ShowerAid™ are trademarks of Atos Medical AB.

For information about protective rights (e.g. patents), please refer to our web page www.atosmedical.com/patents.

Contents

EN - ENGLISH.....	9
DE - DEUTSCH.....	18
FR - FRANÇAIS.....	28
ES - ESPAÑOL.....	38
PT - PORTUGUÊS.....	47
PL - POLSKI.....	56
ORDERING INFORMATION	66

ENGLISH

Intended use

Provox® Coming Home® is an assortment of products and information for newly laryngectomized patients. It provides guidance on product use, lung rehabilitation, and stoma care at home.

Description of Provox® Coming Home®

Provox Coming Home is intended for single patient use.

The products included are selected to be used day and night after a total laryngectomy. See instructions on how to use the different products below. Additional information can be found in Provox Coming Home Booklet.

Contraindications

The devices should only be used in accordance with the Instructions for Use. Patients without the physical, cognitive, or mental ability required to attach, remove or operate the devices themselves, should not use the devices independently and should only use them if they are under sufficient supervision of a clinician or a trained caregiver.

Only valid for HMEs

The devices should not be used by patients with a low tidal volume, as the added dead space (5 ml) may cause CO₂ (Carbon dioxide) retention.

Only valid for Skin Barrier and Adhesive Remover

Hypersensitive or breached skin or known allergies to ingredients.

User assistance information

For additional help or information, please see Contact card provided.

Device lifetime and Disposal

- HMEs are single use devices that are designed to provide pulmonary rehabilitation both day and night, and are to be replaced in line with daily activities.
- The number and type of HMEs needed during a 24 hour period will vary based on the situations encountered and the frequency of involuntary coughing.
- Typically, one HME is used during the night and about two HMEs during the day. The HME may need to be replaced more often if clogged due to mucus secretions or other contamination. HMEs are not reusable and cannot be rinsed in water or any other solution - this washes out the special salt that is essential for HME effectiveness.
- Always follow medical practice and national requirements regarding biohazards when disposing of a used medical device.

Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

Provox® XtraHME™

Provox XtraMoist and Provox XtraFlow

Intended use

Provox XtraHME Cassette is a single use, specialized device intended for patients breathing through a tracheostoma. It is a heat and moisture exchanger (HME) that heats and humidifies inhaled air by retaining heat and moisture from exhaled air in the device. It partially restores lost breathing resistance. For patients with a voice prosthesis or surgical speech fistula it may also facilitate voicing.

Description of the device

The Provox XtraHME Cassettes are single use devices for pulmonary rehabilitation. They are a part of the Provox HME System, which consists of HME cassettes, attachment devices and accessories.

The Provox XtraHME Cassettes have a calcium chloride treated foam in a plastic housing. The top lid can be pushed down with a finger to close the cassette and redirect air through the voice prosthesis to enable speech. After releasing the finger, the top lid returns to its rest position.

The Provox XtraHME Cassettes are available in two versions:

- Provox XtraMoist HME is intended for use during normal everyday activity.
- Provox XtraFlow HME is intended for use during physical activities since it has a lower breathing resistance. It can also be used in a two-step approach to get adapted to the higher breathing resistance of the Provox XtraMoist HME.

WARNINGS

- Be careful not to exert pressure on the lid of the HME unintentionally. Unintentional closing of the lid may cause difficulty in breathing.
- Always inform the patient, caretakers and others about the closing feature of the HME cassette to ensure that they understand its function. Closing the airway to allow voicing is a well-known feature for the laryngectomized patient with a voice prosthesis. For patients without a voice prosthesis or tracheostomized patients this feature might be unknown.

PRECAUTIONS

- Always test the function of Provox XtraHME prior to use. The top lid should immediately return to its open position after releasing the pressure.
- Do not disassemble Provox XtraHME since this will interfere with its proper function.
- Do not reuse Provox XtraHME or attempt to rinse it with water or any other substance. This will substantially reduce the function of the HME. It also increases the risk of potential infections since e.g. bacteria may start to grow in the foam.
- Do not use Provox XtraHME longer than 24 hours. This will increase the risk of potential infections since e.g. bacteria may start to grow in the foam.
- Do not administer medicated nebulizer treatment over the device since the medication can become deposited in the device.
- Do not use humidifiers or heated humidified oxygen via a mask over the tracheostoma while using the device. The HME will become too wet. If oxygen therapy is required, use only non-heated humidified oxygen.

How to use Provox XtraHME

Attach and remove

Place the HME in the holder of your attachment device (e.g. Provox Adhesive or Provox LaryTube) and breathe normally (Fig. 1.1).

If you have a voice prosthesis and you have been cleared to use it for speaking, the lid of Provox XtraHME can be pressed down with your finger to occlude the stoma for speaking. After releasing the finger pressure, the lid opens and you can breathe (Fig. 1.2). Note: Always release the lid completely at inhalation to avoid increased breathing resistance.

To remove the HME, hold the attachment device in place with two fingers and remove the HME from the holder (Fig. 1.3).

Provox® Micron HME™

Intended use

Provox Micron HME is a heat and moisture exchanger (HME) and air filtration device for patients breathing through a tracheostoma. Provox Micron HME partially restores lost breathing resistance. For patients with a voice prosthesis or surgical speech fistula it may also facilitate voicing.

Provox Micron HME is intended to be used with the attachment devices in Provox HME System.

Description of the device

Provox Micron HME is a Heat and Moisture Exchanger combined with an electrostatic filter. The HME is a foam which contains a salt (Calcium Chloride). The HME retains the heat and moisture of the exhaled air. When inhaling, the retained heat and moisture in the HME is given back to the lungs. HME use may help improve the function of the lungs and reduce problems with e.g. coughing and mucus production. New users may experience slight discomfort in the beginning, related to increased breathing resistance. During the first weeks of use, mucus production may seem to increase. This is normal and means that the mucus is getting thinner and easier to cough up.

After a few weeks of HME use, this should stabilize and coughing and mucus production usually decreases.

The Provox Micron HME lid can be pressed down to occlude the stoma in order to speak with a voice prosthesis. When the pressure is released, the lid automatically comes up and the airway passage opens.

Provox Micron HME helps to filter inhaled and exhaled air through consistent normal use. Thereby, small airborne particles, e.g., bacteria, viruses, dust and pollen are restricted from passing through the device into the lungs upon inspiration (see technical data below), while also protecting others nearby during expiration.

Note: Provox Micron is not intended to be used as a Personal Protective Equipment during work that requires breathing protection.

WARNINGS

- Unintentional or accidental pressure on the Provox Micron HME lid may cause difficulty in breathing or suffocation.
- Inform the patient, caregivers and others of this feature to ensure that they understand the closing function. By putting pressure on the Provox Micron HME lid, the airway passage is closed. Closing the airway in order to allow voicing is a well-known feature for the laryngectomized patient with a voice prosthesis, but may be unknown for patients without a voice prosthesis.

PRECAUTIONS

- Provox Micron HME provides good protection through consistent normal-use, as long as there is no air leakage. However, since there are other pathways for e.g. viruses and bacteria to enter the human body, total protection can never be guaranteed.
- Replace Provox Micron HME when needed, the same device must not be used for more than 24 hours after initial use. This can increase the risk for infection due to growth of e.g. bacteria.

- Do not wash and re-use the device. Washing the HME impairs the filtering- and HME functions.
- Do not disassemble Provox Micron HME. Disassembly will destroy its function.
- Do not administer medicated nebulizer treatment over the device since the medication can be deposited in the device.
- Do not use humidifiers or heated humidified oxygen over the device since the HME will become too wet.

How to use Provox Micron HME

Attach and remove

Place Micron HME in the holder of your attachment device (e.g. Provox Adhesive or Provox LaryTube) and breathe normally (Fig. 2.1).

If you have a voice prosthesis and you have been cleared to use it for speaking, the lid of Provox Micron HME can be pressed down with your finger to occlude the stoma for speaking. This will direct the exhaled air through the voice prosthesis. When the pressure is released the airflow will pass through Provox Micron HME again (Fig. 2.2).

To remove Micron HME, hold the attachment device in place with two fingers and remove the HME from the holder (Fig. 2.3).

Provox Micron HME can be removed during or after coughing if the stoma needs to be cleaned from mucus. If you experience any problems with the product, contact your clinic.

Provox® Adhesives

Provox® OptiDerm™, Provox® FlexiDerm™ and Provox® StabiliBase™

Intended use

Provox Adhesives are single use devices intended for laryngectomized patients breathing through a tracheostoma. The devices are attached to the skin around the tracheostoma in order to provide attachment of components of Provox HME System.

Description of the device

Provox Adhesives are adhesives for everyday use, they are designed to ensure an airtight attachment for the Provox HME system components. The adhesives consist of an adapter and an adhesive part with peel-off liners.

Provox FlexiDerm is a very flexible material and has the strongest adhesive properties. It is a sticky, yet soft and flexible adhesive.

Provox OptiDerm is made of a hydrocolloid material that is suitable for sensitive skin. For example, dependent on clinical suitability and continued tolerance, it can be used postoperatively and during and after radiotherapy. The hydrocolloid material is commonly used in wound care and one of its main characteristics is that it absorbs wound fluids. When the material has absorbed fluids it becomes white and depending on the amount of secretions it can become jelly-like and it will no longer stick to the skin. Please note that the material may also absorb mucus coming from the stoma, and therefore it remains very important to always wipe and clear the mucus away from your adhesive and stoma area after you cough.

WARNINGS

- Do not use Provox FlexiDerm, Provox StabiliBase or Provox XtraBase directly after laryngectomy.
- The use of Provox Adhesives need to be assessed both on an individual basis and reviewed regularly through the course of radiotherapy treatment.
- Reuse will affect performance and may cause transfer of micro-organisms leading to infections.
- Only use Provox Adhesives with compatible HMEs and accessories.

PRECAUTIONS

- The adhesive may irritate the skin. Stop using the adhesive if skin irritation develops and consult your clinician.
- When cleaning the skin from residual glue, prevent particles/fluids from entering the tracheostoma.
- When using OptiDerm in the postoperative period or on sensitive skin, it should be removed very slowly and carefully.

How to use a Provox Adhesive

Preparation

Inspect that the adhesive is not damaged or has tears in the adhesive or around the base plate. Always clean the skin with a Provox Cleaning Towel and/or soap and water and let the skin dry before applying the adhesive.

Application

If needed, you can apply skin protection products, like Provox Skin Barrier or adhesive supporting products such as Provox Silicone Glue (read the Instructions for Use accompanying the product). For guidance on how to apply, see figures 3.1 – 3.8 and 4.1 – 4.4.

Removal

Peel the adhesive off gently, using the finger lift tab. An adhesive remover, e.g. Provox Adhesive Remover (read the Instructions for Use accompanying the product), may be helpful for removing adhesives. Always clean the skin with a Provox Cleaning Towel and/or soap and water after removal. Dry the area carefully.

Provox® Skin Barrier

Intended use

Provox Skin Barrier is a single use wipe for laryngectomized patients that forms a barrier between Provox Adhesive and the skin.

WARNINGS

- For external use only.
- Keep out of reach of children.
- Avoid contact with eyes. In the case of accidental contact, flush eyes well with water.
- Do not apply directly to open wounds and mucous membranes.
- Single patient use only. Re-use may cause cross contamination.
- Do not cut or cut the tissue in pieces, small pieces may be aspirated.

PRECAUTIONS

- Avoid deep inhalation of solvent during application and be careful that liquid does not drip into stoma.
- In case you experience signs of irritation or redness, discontinue use of product. If the irritation persists, consult your physician.

How to use a Provox Skin Barrier

Application

Clean the skin around the stoma and let dry. Apply the Provox Skin Barrier to the skin around the stoma that will be covered by the adhesive. Allow the barrier to dry. Apply a Provox Adhesive.

Provox Adhesive Remover

Intended use

Provox Adhesive Remover is a single use wipe to help laryngectomized patients remove Provox Adhesives and Provox Silicone Glue.

WARNINGS

- For external use only.
- Keep out of reach of children.
- Avoid contact with eyes. In the case of accidental contact, flush eyes well with water.
- Do not apply directly to open wounds and mucous membranes.
- Single patient use only. Re-use may cause cross contamination.
- Do not cut the tissue in pieces, small pieces may be aspirated.

PRECAUTIONS

- Avoid deep inhalation of solvent during application and be careful that liquid does not drip into stoma.
- In case you experience signs of irritation or redness, discontinue use of product. If the irritation persists, consult your physician.

How to use Adhesive Remover

Application

Stop the additional oxygen supply, if applied, and remove the HME before applying Adhesive Remover. When using a Provox FlexiDerm, Provox XtraBase, or Provox StabiliBase, first apply the Provox Adhesive Remover on top of the adhesive. This is not necessary for the other adhesives. Then, for all Provox adhesives, grasp the finger lift tab and apply Provox Adhesive Remover on the edge and underneath the adhesive. Carefully remove the adhesive and adhesive residue from the skin. Clean the skin with water and soap or use Provox Cleaning Towel afterwards.

Provox® ShowerAid

Intended use

Provox ShowerAid is used to temporarily replace the HME during showering. The ShowerAid can be placed in all Provox appliance holders.

WARNINGS

- Single patient use only. Reuse between patients may cause cross-contamination.
- Do not bathe or swim, the device does not prevent water from entering the stoma.

How to use a Provox ShowerAid

Attach and remove

Prior to entering the shower, remove the HME and insert Provox ShowerAid with the opening facing down. After showering, remove Provox ShowerAid and insert an HME (Fig. 5.1 – 5.4).

Cleaning and disinfection

Clean the device after each use (Fig. 5.5 – 5.7). Disinfect monthly using 70% Ethanol or Isopropylalcohol for 10 minutes or 3% Hydrogenperoxide for 60 minutes (Fig. 5.8 – 5.9). When the device shows signs of damage, it must be discarded. Replace Provox ShowerAid at least yearly.

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