TRACOE®

technic

Instructions for Use TRACOE® smart Cuff Manager

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ONLY

Federal (USA) law restricts this device to the sale by or on the order of a physician.

TRACOE[®]
technic

USA

Instructions for Use TRACOE[®] smart Cuff Manager REF 730, REF 730-5

NOTE: Please read the Instructions for Use carefully before using the product. This information is considered part of the product and must be available at all times for the individuals responsible for adtstering patient care. When using tubes and accessories the relevant instructions for use must be followed.

CAUTION: The product may only be used by trained physicians and nursing personnel.

1. General Information

1.1 Indications for Use/Intended Use

The TRACOE smart Cuff Manager is intended to maintain the HVLP (high volume low pressure) cuff pressure, of an endotracheal tube or tracheostomy tube, within the range of 20 to 30 cmH₂O through passive control.

It is for single patient use, and indicated for mechanically ventilated or spontaneously breathing patients. It is used under medical supervision in hospitals, pre-hospital (EMS), extended care facilities or outpatient clinics and is suitable for inter- or intra-facility transport.

1.2 Functional Description

The TRACOE smart Cuff Manager includes a spherical, hard shell (transparent) outer casing with an elastic inner balloon (blue), which is attached to a base housing with both a male and female luer connector (figure 1). This device is designed to ensure that the endotracheal/tracheostomy tube cuff pressure remains within the range of 20 cmH₂O to 30 cmH₂O while in clinical use. This range is achieved through passive control, with the inner balloon (blue) of the device acting as a pressure reservoir allowing the flow of air between the balloon and the endotracheal/tracheostomy cuff to maintain the pressure range.

When the endotracheal/tracheostomy tube cuff pressure increases due to compression through patient motion, coughing etc., the surplus air flows slowly into the inner balloon of the TRACOE smart Cuff Manager through the integrated damping function. The damping function reduces sudden changes in cuff pressure which can result in gaps (air/ fluid leakage) between the cuff and the trachea. When the tube cuff pressure is reduced (following the patient motion) the inner balloon will allow the stored air, within the balloon, to freely flow (without damping) back into the endotracheal/tracheostomy tube cuff in order to reinstate the cuff pressure to the prescribed range of 20 cmH₂O to 30 cmH₂O.

CAUTION: The damping function of the inner balloon is only effective in one direction to allow inflation or monitoring of the endotracheal/tracheostomy cuff pressure and therefore, it is not active when a cuff pressure monitor or syringe is inserted in the luer -connector (on the inflation side of the device (figure 1).

CAUTION: The damping function of the TRACOE smart Cuff Manager may result in brief (seconds) increases in cuff pressure before returning to the target range of 20 cmH₂O to 30 cmH₂O.

CAUTION: The endotracheal/tracheostomy tube cuff pressure cannot be maintained if:

- The TRACOE smart Cuff Manager is not initially inflated to a pressure within the range of 20 cmH₂O to 30 cmH₂O and the inner balloon is not expanded to approximately 2/3 to 3/4 of its maximum volume (figure 3).
- There is an air leak in the endotracheal/tracheostomy cuff, inflation line or TRACOE smart Cuff Manager.

CAUTION: The TRACOE smart Cuff Manager must only be used with tubes with HVLP cuffs and a prescribed pressure range of 20 cmH₂O to 30 cmH₂O.

CAUTION: The cuff inflation line of the endotracheal/tracheostomy tube can be damaged by sudden and powerful pulling on the tube therefore the TRACOE smart Cuff Manager should be placed such that it maintains free movement in order to avoid unintentional snagging or pulling.

2. General Description

The TRACOE smart Cuff Manager is designed to ensure that the endotracheal/tracheostomy tube HVLP cuff pressure remains within the range of 20 cmH₂O to 30 cmH₂O while in clinical use. This

range is achieved through passive control, with the inner balloon (blue) of the device acting as a pressure reservoir allowing the flow of air between the balloon and the endotracheal/tracheostomy cuff to maintain the pressure range.

The elastic inner balloon (blue) is located inside a spherical hard shell (transparent) outer casing. This inner balloon is connected to the base housing (white), which has two luer-connectors (figure 1):

- A male luer-connector on the cuff side serves to connect the TRACOE smart Cuff Manager to the cuff inflation line (pilot balloon) of an endotracheal/tracheostomy tube.
- A female luer-connector on the opposite (inflation) side serves to connect the TRACOE cuff pressure monitor (cpm), or any other suitable manometer/ inflation device (syringe).


Fig. 1: Components and connectors of the TRACOE smart Cuff Manager

The female luer-connector (inflation side) is used for:

- inflating / re-inflating the inner balloon of the TRACOE smart Cuff Manager,
- filling the endotracheal/tracheostomy tube cuff,
- measuring the cuff pressure at regular intervals or
- re-inflating the inner balloon (as required).

This female luer connector has an incorporated self-sealing valve which prevents air flow upon removal of a syringe or pressure measuring unit, see figure 1.

Packaging: The TRACOE smart Cuff Manager is provided sterile (sterilized with ethylene oxide). This allows the product to be connected under sterile conditions. Sterility is only ensured if the sterile device packaging is intact and sealed on all sides prior to opening.

Material composition: ABS, PETG, TPE, silicone and stainless steel. The product is not made with natural rubber latex or phthalates (e.g. DEHP).

Disposal: Following clinical use of the TRACOE smart Cuff Manager, the product should be considered a bio-hazardous material and disposed of in accordance with hospital and/or national regulations.

Dimensions and weight:

Length: 55 mm, Height: 77 mm, Weight: 16 g

Storage: The product must be stored in a cool, dry place, protected from sunlight and extreme cold (10°C to 30°C).

Maximum Useful Life: The maximum period of clinical use is 29 days. The device must be replaced if it is damaged or does not meet specification.

Product for Single-Patient Use: This device can be connected and reconnected multiple times on a single patient over the maximum period of clinical use.

Transport: The TRACOE smart Cuff Manager can remain connected during inter- or intra-facility transport.

Environment: Normal, weather-related fluctuations in air pressure do not need to be considered when using the TRACOE smart Cuff Manager.

Limitations: The TRACOE smart Cuff Manager does not control the cuff pressure to a specific set point but through passive control it maintains the cuff pressure within the established safety limit of 20 cmH₂O to 30 cmH₂O.

Packaging Contents:

- REF 730 1 pc. TRACOE smart Cuff Manager (individually packed, sterile)
- REF 730-5 5 pcs TRACOE smart Cuff Manager (individually packed, sterile)

3. Indications

TRACOE smart Cuff Manager is suitable for patients where sealing of the trachea by means of an HVLP cuffed endotracheal/tracheostomy tube (with an indicated cuff pressure range of 20 cmH₂O to 30 cmH₂O) is required. Optimal cuff pressure for each patient can only be set by the attending caregiver.

It is for single patient use, and indicated for mechanically ventilated or spontaneously breathing patients. It is used under medical supervision in hospitals, pre-hospital (EMS), extended care facilities or outpatient clinics and is suitable for inter- or intra-facility transport.

4. Contraindications

- TRACOE smart Cuff Manager is not suitable if the prescribed endotracheal/tracheostomy cuff pressure is greater than 30 cmH₂O, or less than 20 cmH₂O.
- TRACOE smart Cuff Manager is not suitable for use with laryngeal masks.
- TRACOE smart Cuff Manager is considered MR Unsafe and must be removed before entering a magnetic resonance imaging (MRI) area due to a metal spring incorporated in the device.

5. General Precautions

- Prevent fluids from entering the TRACOE smart Cuff Manager ventilation opening (top of the outer casing, figure 1); e.g. do not bath or shower with the device attached.
- Ensure water/fluids do not enter the endotracheal/tracheostomy cuff or inflation line as this can impair the damping function and air flow required to maintain cuff pressure.
- Regular monitoring of the endotracheal/tracheostomy cuff pressure (e.g. every 8 hours) is required as directed by hospital policy or national guidelines.

6. Warnings

- The TRACOE smart Cuff Manager ventilation opening (top of the outer casing) must remain open and unobstructed at all times (figure 1).
- Ensure the endotracheal/tracheostomy inflation filling tube (pilot balloon) is not kinked to allow air flow and maintain cuff pressure between the TRACOE smart Cuff Manager and the endotracheal/tracheostomy tube cuff.

7. Adverse Side Effects

There are no known adverse side effects with this device.

8. Instructions for Use

8.1 Connection and Inflation

With the endotracheal/tracheostomy tube (deflated or inflated HVLP cuff) positioned within the patient, the TRACOE smart Cuff Manager is connected to the

open end of the cuff inflation line (pilot balloon) by inserting the male Luer-connector (Luer m, figure 1), at the base of the TRACOE smart Cuff Manager, and gently turning to secure the connection (figure 2).


Fig. 2: Connect the inflation line to the male Luer-connector

Once the TRACOE smart Cuff Manager is connected to the cuff inflation line, insert the cuff pressure monitor or syringe into the female Luer connection, at the base (Luer f, figure 1), and inject sufficient air to properly inflate the inner blue balloon to a diameter of 2/3 – 3/4 (figure 3) as compared to the outer housing.


Fig. 3: Required inflation size

Following the initial inflation of the TRACOE smart Cuff Manager the device should be left undisturbed for a few minutes (less than 5) in order to equalize the pressure between the TRACOE smart Cuff Manager and the endotracheal/tracheostomy tube cuff.

After this initial rest period, the pressure within the TRACOE smart Cuff Manager must be measured using a cuff pressure monitor. If the pressure, within the TRACOE smart Cuff Manager is above or below the acceptable range of 20 cmH₂O to 30 cmH₂O remove/add enough air to ensure the inner balloon is filled to the 2/3 to 3/4 level and then repeat the rest period before verifying the pressure range again.

Warning

If the TRACOE smart Cuff Manager deflates or does not maintain the 2/3 to 3/4 fill level it may indicate

an air leak in the device, the cuff, or the filling line. It is recommended to replace the TRACOE smart Cuff Manager and if the situation reoccurs, consider replacing the endotracheal/tracheostomy tube.

Warning

Ensure that the endotracheal/tracheostomy tube cuff is properly positioned prior to filling the TRACOE smart Cuff Manager.

When an extension tube (REF 707, REF 702) is required between the cuff inflation line of the endotracheal/tracheostomy tube and the TRACOE smart Cuff Manager, insert the male end of the extension tube into the cuff inflation line, then insert the male Luer connector, of the TRACOE smart Cuff Manager, into the open end of the extension tube and gently turn to secure the connection. Refer to section 9 Accessories for more information on the extension tubes that are approved for use with the TRACOE smart Cuff Manager.

CAUTION: The TRACOE smart Cuff Manager should not be subjected to unnecessary external forces. Do not pull the tubes to disconnect. For more information on disconnecting the TRACOE smart Cuff Manager, section 8.6.

8.2 Inspection of the Inflation Diameter/ Cuff Pressure

During clinical use the TRACOE smart Cuff Manager should be monitored (visually and with a cuff pressure monitor) at regular intervals as defined by clinical procedures (e.g. every 8 hours) to ensure the device maintains the endotracheal/tracheostomy tube cuff pressure. For visual inspection of the device, it was specifically designed with a transparent outer shell and a blue inner balloon which clearly demonstrates the balloon inflation level.

CAUTION: If the blue inner balloon is pressing tightly against or attached to the transparent outer housing it cannot maintain the cuff pressure in the prescribed range. If this occurs, the TRACOE smart Cuff Manager should be disconnected and endotracheal/tracheostomy tube cuff pressure should be checked, using a cuff pressure monitor, to ensure that the cuff pressure is within range. If necessary, refill the endotracheal/tracheostomy tube cuff to the desired value then reconnect the TRACOE smart Cuff Manager and refill the inner balloon to the appropriate diameter, following the procedures described in 8.1.

8.3 Placement of the TRACOE smart Cuff Manager

The TRACOE smart Cuff Manager's lightweight design allows it to be placed directly on or next to the patient during clinical use. If a different location is preferred, extension tubes are available in 0.3 and 1.0 meter lengths (REF 707, REF 702) which can be connected between the cuff inflation line of the endotracheal/tracheostomy tube and the TRACOE smart Cuff Manager, see section 9 Accessories.

Warning

The TRACOE smart Cuff Manager ventilation opening (top of the outer casing) must remain open and unobstructed at all times, see figure 1.

Warning

Ensure that the cuff inflation line of the endotracheal/tracheostomy tube and any connected extension tubing is not kinked to allow air flow and maintain endotracheal/tracheostomy cuff pressure.

8.4 Refilling

Throughout the 29 days of clinical use the inner balloon of the TRACOE smart Cuff Manager may require refilling of air to maintain the 2/3 – 3/4 fill level (figure 3). To refill the inner balloon, follow the inflation procedure defined in section 8.1. If the volume of the balloon is too low, the cuff pressure may be increased.

Warning

If the TRACOE smart Cuff Manager deflates or does not maintain the 2/3 to 3/4 fill level it may indicate an air leak in the device, the cuff, or the filling line. It is recommended to replace the TRACOE smart Cuff Manager and if the situation reoccurs, consider replacing the endotracheal/tracheostomy tube.

8.5 Deflation of the Endotracheal/ Tracheostomy Tube Cuff

The TRACOE smart Cuff Manager must be disconnected from the endotracheal/tracheostomy tube cuff inflation line (pilot balloon) prior to deflation of the HPLV cuff, see section 8.6.

CAUTION: Do not deflate the HPLV cuff while the TRACOE smart Cuff Manager is connected.

8.6 Disconnecting TRACOE smart Cuff Manager

To disconnect, grasp the TRACOE smart Cuff Manager and gently turn while pulling the device from the endotracheal/tracheostomy tube cuff inflation

line. The inner balloon will deflate once the TRACOE smart Cuff Manager is disconnected.

If the TRACOE smart Cuff Manager is deliberately or accidentally disconnected, the endotracheal/tracheostomy cuff should retain its cuff pressure due to the pilot balloon valve of the cuff inflation line.

Note: If reconnecting the TRACOE smart Cuff Manager (with a deflated inner balloon) the endotracheal/tracheostomy tube cuff will briefly lose pressure until the inner balloon of the TRACOE smart Cuff Manager has been refilled, see section 8.1. In this situation the standard clinical procedure for a deflated endotracheal/tracheostomy cuff should be followed, e.g. secretion suction.

9. Accessories

TRACOE *technic* Connection Tube

REF 707 (0.3 m), REF 702 (1.0 m)

These tubes provide a secure extension between the TRACOE smart Cuff Manager and the endotracheal/tracheostomy cuff inflation line. The extension tubes, REF 707 and REF 702 are provided sterile (sterilized with ethylene oxide).

TRACOE *technic* cpm REF 720

This cuff pressure monitor can be used to inflate the TRACOE smart Cuff Manager and measure/adjust the cuff pressure during clinical use.

10. Cleaning and Disinfection

The TRACOE smart Cuff Manager can be wiped clean using a moist cloth with or without mild cleaning agents. If necessary, the outside housing of the device can be wiped with a disinfectant. No other maintenance is required.

CAUTION: The TRACOE smart Cuff Manager cannot be immersed in water, disinfection solution or any other liquid. Care must be taken to avoid fluids from entering the connecting luers and the ventilation opening, see figure 1, as this will impair the damping function and inhibit air flow required to maintain cuff pressure.

11. Information for Ordering

REF Number	Contents	Device
REF 730	1 pc.	TRACOE smart Cuff Manager
REF 730-5	5 pcs.	TRACOE smart Cuff Manager
REF 707	5 pcs.	Connecting tube 0.3 metres
REF 702	10 pcs.	Connecting tube 1.0 metre

12. Returns

The return of used products will only be accepted if a decontamination certificate is enclosed. This form is available either directly from TRACOE medical GmbH, or www.tracoe.com.

13. General Terms and Conditions

The sale, delivery and return of all TRACOE products shall comply with our General Terms and Conditions which are available from TRACOE medical GmbH or www.tracoe.com.

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Instructions for Use TRACOE® smart Cuff Manager REF 730, REF 730-5

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CAUTION: The product may only be used by individuals trained in the use of tracheostomy and endotracheal tubes.

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It is for single patient use, and indicated for mechanically ventilated or spontaneously breathing patients. It is used in hospitals, pre-hospitals (EMS), extended care facilities, outpatient clinics, ambulant care facilities or homecare, and is suitable for inter- or intra-facility transport.

1.2 Functional Description:

The TRACOE smart Cuff Manager includes a spherical, hard shell (transparent) outer casing with an elastic inner balloon (blue), which is attached to a base housing with both a male and female luer connector (figure 1). This device is designed to ensure that the endotracheal/tracheostomy tube cuff pressure remains within the range of 20 cmH₂O to 30 cmH₂O while in use. This range is achieved through passive control, with the inner balloon (blue) of the device acting as a pressure reservoir allowing the flow of air between the balloon and the endotracheal/tracheostomy cuff to maintain the pressure range.

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2. General Description

The TRACOE smart Cuff Manager is designed to ensure that the endotracheal/tracheostomy tube HVLP cuff pressure remains within the range of 20 cmH₂O to 30 cmH₂O while in use. This range

is achieved through passive control, with the inner balloon (blue) of the device acting as a pressure reservoir allowing the flow of air between the balloon and the endotracheal/tracheostomy cuff to maintain the pressure range.

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Disposal: Following use of the TRACOE smart Cuff Manager, the product should be considered a bio-hazardous material and disposed of in accordance with hospital and/or national regulations.

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Maximum Useful Life: The maximum period of use is 29 days. The device must be replaced if it is damaged or does not meet specification.

Product for Single-Patient Use: This device can be connected and reconnected multiple times on a single patient over the maximum period of use.

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It is for single patient use, and indicated for mechanically ventilated or spontaneously breathing patients. It is used in hospitals, pre-hospitals (EMS), extended care facilities, outpatient clinics, ambulant care facilities or homecare, and is suitable for inter- or intra-facility transport.

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Once the TRACOE smart Cuff Manager is connected to the cuff inflation line, insert the cuff pressure monitor or syringe into the female Luer connection, at the base (Luer f, figure 1) and inject sufficient air to properly inflate the inner blue balloon to a diameter of 2/3 – 3/4 (figure 3) as compared to the outer housing.


Fig. 3: Required inflation size

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During use the TRACOE smart Cuff Manager should be monitored (visually and with a cuff pressure monitor) at regular intervals as defined by clinical procedures (e.g. every 8 hours) to ensure the device maintains the endotracheal/tracheostomy tube cuff pressure. For visual inspection of the device, it was specifically designed with a transparent outer shell and a blue inner balloon which clearly demonstrates the balloon inflation level.

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8.4 Refilling

Throughout the 29 days of use the inner balloon of the TRACOE smart Cuff Manager may require refilling of air to maintain the 2/3 – 3/4 fill level (figure 3). To refill the inner balloon, follow the inflation procedure defined in section 8.1. If the volume of the balloon is too low, the cuff pressure may be increased.

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TRACOE *technic* cpm REF 720

This cuff pressure monitor can be used to inflate the TRACOE smart Cuff Manager and measure/adjust the cuff pressure during use.

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